UNICORN 7.3.0 1 (14) User: AM/c304117 4/2/2025 9:45:13 AM -04:00 Method: v001 Scouting Method LB2273 CEX Text instructions Main method: Base: CV, Vc=69.555 {1}, LB2273 CEX 24pt6cm H x 60cm D 0.00 Phase: Method Settings Base: SameAsMain 0.00 0.00 Phase: User Defined Base: SameAsMain 0.00 Set mark: (Result Name) #Result Name 0.00 0.00 Block: Start Conditions Base: SameAsMain 0.00 0.00 Air Alarm: Disabled, Disabled Comment: See valley detect method if using valley strategy in elution. Comment: Disabling Flow warning prevents "no flow" error when method is run. 0.00 Flow warning: Disabled Comment: Reduce flow deviation if method will require flowrate less than 10 L/hr FlowDeviation FIT PA: 20.0 {1/hour}, -20.0 {1/hour}, 300.0 {sec}, Enabled FlowDeviation FIT PB: 20.0 {1/hour}, -20.0 {1/hour}, 300.0 {sec}, Enabled 0.00 PIT PA: 5.00 {bar}, 0.00 {bar}, 3.50 {bar}, 0.00 {bar}, 0.00 {bar}, Enabled 0.00 PIT PB: 3.00 {bar}, 0.00 {bar}, 2.70 {bar}, 0.00 {bar}, 0.00 {bar}, Enabled 0.00 Wavelength: 280 {nm}, 0 {nm}, 0 {nm} 0.00 End Block Comment: THROUGHOUT: Update inlet purges to 7L, 10L, and 15L for 3/8", 1/2", and 3/4" skid respectively Comment: THROUGHOUT: Manflow = 60% for 3/4" skid. Manflow = 100% for 3/8" and 1/2" skid Comment: THROUGHOUT: Ensure totalizer reset matches pump in block Block: (Startup Blocks) #Startup Blocks 0.00 Base: SameAsMain 0.00 Block: Prepare Purge Col Bypass 0.00 Base: Time, ColumnSameAsMain 0.00 Message: Prepare to purge inlet 5 and air from bypass line and installed hoses. Ensure column is bypassed. , Screen, No sound 0.00 Pause: Infinite {min} 0.01 End Block 0.00 Block: Purge Inlet 5 Col Bypass 0.00 Base: Volume, ColumnSameAsMain 0.00 Air Alarm: Disabled, Disabled 0.00 Inlet: Closed, Inlet5 0.00 BubbleTrap: Bypass

```
UNICORN 7.3.0
                                                                      2 (14)
User: AM/c304117 4/2/2025 9:45:13 AM -04:00
Method: v001 Scouting Method LB2273 CEX
        0.00
               Filter: Bypass
        0.00
               Column: UpFlow
        0.00
               Outlet: Waste
        0.00
               ManFlow: 60.0 {%}
        0.00
               Set mark: Purge Inlet 5 and Column Bypass Line: QD00008
        Comment: Breakpoint settings: 10L up and down for 3/8" and 1/2", 15L
        up and down for 3/4"
        15.00
                Column: DownFlow
        30.00
                End Block
      0.00
             Block: Place Column Inline
               Base: Time, ColumnSameAsMain
        0.00
               Message: Close valves to column bypass loop and open column
        0.00
        valves (column inline)., Screen, No sound
        0.00
               Pause: Infinite {min}
        0.01
               End Block
      0.00
             Block: Purge Inlet Sample
        0.00
               Base: Volume, ColumnSameAsMain
               Air Alarm: Disabled, Disabled
        0.00
        0.00
             Inlet: Sample, Closed
        0.00
               BubbleTrap: Bypass
        0.00
              Filter: Bypass
               Column: Bypass Both
        0.00
        0.00
             Outlet: Waste
        0.00
             ManFlow: 60.0 {%}
        0.00
               Set mark: Purge Inlet Sample: QD00305
        15.00
               End Block
      0.00
             Block: Connect Equil to Inlet 1
               Base: Time, ColumnSameAsMain
        0.00
        0.00
               Message: Connect QD00305 to Inlet 1 and open clamps., Screen,
        No sound
        0.00
               Pause: Infinite {min}
        0.01
               End Block
      0.00
            Block: Purge Inlet 4
               Base: Volume, ColumnSameAsMain
        0.00
        0.00
               Air Alarm: Disabled, Disabled
        0.00
               Inlet: Closed, Inlet4
        0.00
              BubbleTrap: Bypass
        0.00
               Filter: Bypass
        0.00
               Column: Bypass Both
        0.00
               Outlet: Waste
        0.00
              ManFlow: 60.0 {%}
        0.00
               Set mark: Purge Inlet 4: QD00007
               End Block
        15.00
      0.00
             Block: Purge Inlet 7
```

```
UNICORN 7.3.0
                                                                      3 (14)
User: AM/c304117 4/2/2025 9:45:13 AM -04:00
Method: v001 Scouting Method LB2273 CEX
        0.00
               Base: Volume, ColumnSameAsMain
               Air Alarm: Disabled, Disabled
        0.00
        0.00
               Inlet: Closed, Inlet7
        0.00
               BubbleTrap: Bypass
               Filter: Bypass
        0.00
        0.00
               Column: Bypass Both
        0.00
               Outlet: Waste
        0.00
             ManFlow: 60.0 {%}
        0.00
               Set mark: Purge Inlet 7: QD00346
        15.00
               End Block
      0.00
             End Block
         Block: (Prepare Purge Skid Filter) #Prepare Purge Skid Filter
             Base: Time, ColumnSameAsMain
      0.00
             Message: Prepare to purge the skid filter., Screen, No sound
      0.00
             Pause: Infinite {min}
      0.01
             End Block
    0.00
          Block: (Purge Skid Filter Inlet 1) #Purge Skid Filter Inlet 1
             Base: Volume, ColumnSameAsMain
      0.00
             Air Alarm: Disabled, Disabled
      0.00
      0.00 Inlet: Inlet1, Closed
      0.00 BubbleTrap: Inline
      0.00
           Filter: Inline
      0.00 Column: Bypass Both
      0.00 Outlet: Waste
      0.00
           ManFlow: 60.0 {%}
      0.00
             Set mark: Purge Inlet 1 and Filter
      Comment: Regardless of skid size, keep 20L breakpoint with standard
      skid filter setup.
      Comment: Adjust breakpoint in scenarios with non-standard setups (i.e.
      glass fiber + cartridge filter).
      20.00
              End Block
    0.00
           Block: (Flush Outlet Mainstreams Equil)
    #Flush Outlet Mainstreams Equil
             Base: Volume, ColumnSameAsMain
      0.00
      0.00
             Air Alarm: Disabled, Disabled
      0.00
             Inlet: Inlet1, Closed
      0.00
           BubbleTrap: Inline
      0.00
           Filter: Inline
      0.00
           Column: Bypass Both
      0.00
           ManFlow: 60.0 {%}
      0.00
            Fractions: 3, 5.0 {1}, Outlet1
      Comment: Set breakpoint to number of mainstreams times 5L.
      Comment: Always keep with filter flush; block strategically placed to
      provide additional filter flush.
```

```
UNICORN 7.3.0
                                                                     4 (14)
User: AM/c304117 4/2/2025 9:45:13 AM -04:00
Method: v001 Scouting Method LB2273 CEX
      15.00
              End Block
    0.00
         Block: (Pause attach outlet containers)
    #Pause attach outlet containers
      0.00
             Base: Time, ColumnSameAsMain
             Message: Attach Outlets to effluent containers per ticket
      0.00
      instructions., Screen, No sound
           Pause: Infinite {min}
      0.01
             End Block
         Block: (Pre-Use Sanitization) #Pre-Use Sanitization
             Base: SameAsMain
      0.00
             Block: Pre Use Sani
        0.00
               Base: SameAsMain
        0.00
               Air Alarm: Disabled, Enabled
        0.00
             FIT PB Totalizer Reset
        0.00 Inlet: Closed, Inlet4
        0.00
               BubbleTrap: Inline
        0.00
             Filter: Inline
        0.00
               Column: DownFlow
        0.00
             Outlet: Waste
        0.00 Flow: (300) #Pre Sani Flowrate {cm/h}
        0.00 Set mark: Pre Sani
        2.00
               Snapshot: Pre Sani End
        2.00
               End Block
      0.00
           Block: Pause Complete Sanitization
        0.00
             Base: Time, ColumnSameAsMain
        0.00
               Pause: 30.00 {min}
        0.01
               End Block
             End Block
      0.00
    0.00 Block: Purge A Pump
      0.00 Base: Time, ColumnSameAsMain
      0.00
             Air Alarm: Disabled, Disabled
      0.00
             Inlet: Inlet1, Closed
      Comment: Change inlet to match first inlet used from pump A,
      henceforth (from this point forward)
      Comment: Inlet purge volumes do not apply to pump purges and block
      should remain in base time.
      0.00
            BubbleTrap: Bypass
      0.00 Filter: Bypass
      0.00 Column: Bypass Both
           Outlet: Waste
      0.00
      0.00 ManFlow: 60.0 {%}
      0.00
             Set mark: Purge Inlet 1
      2.00
             End Block
           Block: Purge B Pump
```

```
UNICORN 7.3.0
                                                                     5 (14)
User: AM/c304117 4/2/2025 9:45:13 AM -04:00
Method: v001 Scouting Method LB2273 CEX
             Base: Time, ColumnSameAsMain
      0.00
      0.00
             Air Alarm: Disabled, Disabled
             Inlet: Closed, Inlet7
      0.00
      Comment: Change inlet to match first inlet used from pump B,
      henceforth (from this point forward)
      Comment: Inlet purge volumes do not apply to pump purges and block
      should remain in base time.
      0.00
             BubbleTrap: Inline
      0.00 Filter: Inline
      0.00 Column: Bypass Both
      0.00 Outlet: Waste
      0.00 ManFlow: 60.0 {%}
      0.00 Set mark: Purge Inlet 7
      2.00
             End Block
    0.00
         Block: Pre Equilibration
      0.00
             Base: SameAsMain
      0.00
             Air Alarm: Disabled, Enabled
      0.00 FIT PB Totalizer Reset
      0.00 Inlet: Closed, Inlet7
      0.00 BubbleTrap: Inline
      0.00 Filter: Inline
      0.00
           Column: DownFlow
      0.00 Outlet: Waste
      0.00 Flow: (300) #Equil Flowrate {cm/h}
      0.00 Set mark: Pre Equil
      2.00
             Snapshot: Pre Equil End
      2.00
             End Block
    0.00
         Block: Equilibration
      0.00
             Base: SameAsMain
      0.00
             Air Alarm: Disabled, Enabled
      0.00
             FIT_PA_Totalizer_Reset
      0.00
             FIT PB Totalizer Reset
            FIT PA+PB Totalizer Reset
      0.00
      0.00 Inlet: Inlet1, Closed
      0.00
            BubbleTrap: Inline
      0.00 Filter: Inline
      0.00
            Column: DownFlow
      0.00 Outlet: Waste
      0.00
           Flow: (300) #Equil Flowrate {cm/h}
      0.00 Set mark: Equil
      3.00
             Snapshot: Equil End
      3.00
             End Block
    0.00
          Block: UV Auto Zero
             Base: Time, ColumnSameAsMain
```

```
UNICORN 7.3.0
                                                                     6 (14)
User: AM/c304117 4/2/2025 9:45:13 AM -04:00
Method: v001 Scouting Method LB2273 CEX
             AT PF AZ
      0.00
             End Block
      0.10
           Block: (Connect_Charge_to_Inlet_Sample)
    #Connect Charge to Inlet Sample
            Base: Time, ColumnSameAsMain
      0.00
             Message: Connect Charge to Inlet Sample and open clamps.,
      Screen, No sound
      0.00
             Pause: Infinite {min}
      0.01
             End Block
    0.00
         Block: Charge
      0.00
           Base: Volume, ColumnSameAsMain
      0.00
             Injection Mark
      0.00 Air Alarm: Disabled, Enabled
      0.00 FIT PA Totalizer Reset
      0.00 Inlet: Sample, Closed
      0.00 BubbleTrap: Inline
      0.00 Filter: Inline
           Column: DownFlow
      0.00
      0.00
           Outlet: Waste
      0.00
           Flow: (200) #Charge Flowrate {cm/h}
      0.00 Set mark: Charge
      (800.00) #Set Charge Volume Snapshot: Charge End
      800.00
             End Block
    0.00
           Block: Post Charge Wash
             Base: SameAsMain
      Comment: Ensure first column wash flows through the material load
      pump.
      0.00
             Air Alarm: Disabled, Enabled
      0.00 FIT PA Totalizer Reset
      0.00 FIT PB Totalizer Reset
      0.00 FIT PA+PB Totalizer Reset
      0.00
             Inlet: Inlet1, Inlet7
      0.00
           BubbleTrap: Inline
           Filter: Inline
      0.00
      0.00
            Column: DownFlow
           Outlet: Waste
      0.00
      0.00
           GradMode: FlowGradient
      0.00 Gradient: 20.0 {%B}, 0.00 {base}
      0.00
             Flow: (200) #Post Charge Wash Flowrate {cm/h}
             Set mark: Post Charge Wash
      0.00
      2.00
             Snapshot: Post Charge Wash End
      2.00
             Watch off: AT PF UV 1
      2.00
             End Block
    0.00
           Block: (Elution) #Elution
```

```
UNICORN 7.3.0
                                                                      7 (14)
User: AM/c304117 4/2/2025 9:45:13 AM -04:00
Method: v001 Scouting Method LB2273 CEX
      0.00
             Base: SameAsMain
      0.00
           Air Alarm: Disabled, Enabled
           FIT PA Totalizer_Reset
      0.00
      0.00 FIT PB Totalizer Reset
      0.00 FIT PA+PB Totalizer Reset
      0.00 Inlet: Inlet1, Inlet7
      0.00
           BubbleTrap: Inline
      0.00 Filter: Inline
           Column: DownFlow
      0.00
      0.00 Outlet: Waste
      0.00
           GradMode: FlowGradient
      0.00 Flow: (200) #Elution Flowrate {cm/h}
      0.00 Gradient: 70.0 {%B}, 12.50 {base}
      0.00
           Set mark: Elution
      0.00
           Watch: PIT PA., Greater than, 3.00 {bar}, Elution Message
             Base: SameAsMain
        0.00
             Message: Pause during elution may trigger UV and impact
        0.00
        mainstream collection, Screen, No sound
        0.00
               End Block
            Watch: PIT PB., Greater than, 2.70 {bar}, Slow Flow
      0.00
               Base: SameAsMain
               Flow: (150) #Reduced Elution Flow {cm/h}
        Comment: Ensure Reduced Elution Flowrate is adequately below defined
        flowrate.
               Message: Pause during elution may trigger UV and impact
        mainstream collection, Screen, No sound
               End Block
        0.00
      Comment: LHM4350 2mm Path Length Compensation Factor = 4.92
             Block: Watch UV
        0.00
               Base: SameAsMain
        0.00
               Watch: AT PF UV 1, Greater than, 0.6098 {AU}, Collect Peak
                 Base: SameAsMain
          0.00
          0.00
                 Injection Mark
          0.00
                 Outlet: (Outlet1) #MS Outlet
                 Watch: AT PF UV 1, Greater than, 1.5000 {AU},
          watch Less Than
            0.00
                  Base: SameAsMain
            Comment: Design peak protect to be adequately above backside cut
            but within skid capability (ex. ~NMT 3 AU).
                   Watch: AT PF UV 1, Less than, 0.8130 {AU}, End block
            0.00
                   End Block
            0.00
          Comment: BS cut triggers end block command, ends elution block, &
          automatically transitions to next block.
          0.00
                 End Block
        0.00
               End Block
```

```
UNICORN 7.3.0
                                                                     8 (14)
User: AM/c304117 4/2/2025 9:45:13 AM -04:00
Method: v001 Scouting Method LB2273 CEX
      12.50
              GradMode: FlowGradient
      12.50
              Gradient: 70.0 {%B}, 0.00 {base}
            Snapshot: Elution End
      14.50
      14.50 Watch off: AT PF UV 1
      14.50
            End Block
      Comment: Keep following gradient/gradmode commands if PFC has
      isocratic hold and update per PFC requirements
           Block: (Elution Full Backside) #Elution Full Backside
    0.00
      0.00
             Base: SameAsMain
      0.00
             Air Alarm: Disabled, Enabled
      0.00
           FIT PA Totalizer Reset
      0.00 FIT PB Totalizer Reset
      0.00 FIT PA+PB Totalizer Reset
      0.00 Inlet: Inlet1, Inlet7
      0.00 BubbleTrap: Inline
      0.00 Filter: Inline
      0.00 Column: DownFlow
      0.00 Outlet: Waste
      0.00 Flow: (200) #Elution Flowrate {cm/h}
      0.00 GradMode: FlowGradient
      0.00 Gradient: 70.0 {%B}, 12.50 {base}
      0.00
           Set mark: Elution
      0.00
           Watch: PIT PA., Greater than, 3.00 {bar}, Elution Message
        0.00 Base: SameAsMain
             Message: Pause during elution may trigger UV and impact
        mainstream collection, Screen, No sound
               End Block
        0.00
      Comment: LHM4350 2mm Path Length Compensation Factor = 4.92
             Watch: PIT PB., Greater than, 2.70 {bar}, Slow Flow
        0.00
               Base: SameAsMain
        0.00
               Flow: (150) #Reduced Elution Flow {cm/h}
        Comment: Ensure Reduced Elution Flowrate is adequately below defined
        flowrate.
        0.00
               Message: Pause during elution may trigger UV and impact
        mainstream collection, Screen, No sound
               End Block
      0.00
            Block: Watch UV 1
        0.00
               Base: SameAsMain
               Watch: AT_PF_UV_1, Greater than, 0.6098 {AU}, Collect Peak 1
        0.00
          0.00
                 Base: SameAsMain
          0.00
                 Injection Mark
          0.00
                 Outlet: (Outlet1) #MS Outlet
                 Watch: AT PF UV 1, Greater than, 1.5000 {AU},
          watch Less Than 1
            0.00
                   Base: SameAsMain
```

```
UNICORN 7.3.0
                                                                      9 (14)
User: AM/c304117 4/2/2025 9:45:13 AM -04:00
Method: v001 Scouting Method LB2273 CEX
            Comment: Design peak protect to be adequately above backside cut
            but within skid capability (ex. NMT3AU).
            Comment: Design peak protect to be adequately above backside cut
            but within skid capability (ex. ~NMT 3 AU).
                   Watch: AT PF UV 1, Less than, 0.8130 {AU}, stop Collect
                     Base: SameAsMain
               0.00
               0.00
                     Outlet: Waste
              Comment: stop Collect block ends MS collection and sends
               remaining elution to waste.
               0.00
                    End Block
            0.00 End Block
          0.00
                End Block
        0.00
               End Block
      Comment: Keep following gradient/gradmode commands if PFC has
      isocratic hold and update per PFC requirements
              GradMode: FlowGradient
      12.50
      12.50
              Gradient: 70.0 {%B}, 0.00 {base}
      14.50
            Snapshot: Elution End
      14.50
             Watch off: AT PF UV 1
      14.50
              End Block
    0.00
           Block: High Salt Wash
      0.00
             Base: SameAsMain
      Comment: Delete High Salt Wash block if not needed per PFC
      Comment: Turn off the pressure watches from Elution.
            Watch off: PIT PA.
      0.00
      0.00
             Watch off: PIT PB.
             Watch off: AT PF UV 1
      0.00
             Air Alarm: Disabled, Enabled
      0.00
      0.00 FIT PB Totalizer Reset
      0.00 Inlet: Closed, Inlet7
      0.00 BubbleTrap: Inline
      0.00
           Filter: Inline
      0.00 Column: UpFlow
      0.00 Outlet: Waste
      0.00 Flow: (300) #High Salt Wash Flowrate {cm/h}
      0.00
           Set mark: High Salt Wash
             Snapshot: High Salt Wash End
      2.00
      2.00
             End Block
    0.00
         Block: Post Use Sanitization
             Base: SameAsMain
      0.00
      0.00
             Block: Post Use Sani
               Base: SameAsMain
        Comment: Turn off the pressure watches from Elution.
               Watch off: PIT PA.
        0.00
```

```
UNICORN 7.3.0
                                                                     10 (14)
User: AM/c304117 4/2/2025 9:45:13 AM -04:00
Method: v001 Scouting Method LB2273 CEX
               Watch off: PIT PB.
        0.00
               Watch off: AT PF UV 1
        0.00
               New chromatogram: Post Use Sani
        0.00
        0.00
               FIT PB Totalizer Reset
        0.00
               Inlet: Closed, Inlet4
        0.00
             BubbleTrap: Inline
        0.00
               Filter: Inline
        0.00
             Column: UpFlow
        0.00
             Outlet: Waste
        0.00 Flow: (300) #Post Use Sani Flowrate {cm/h}
        0.00
               Set mark: Post Use Sani
        2.00 FlowDeviation FIT PA: 20.0 {1/hour}, -20.0 {1/hour}, 300.0
        {sec}, Disabled
              FlowDeviation FIT PB: 20.0 {1/hour}, -20.0 {1/hour}, 300.0
        2.00
        {sec}, Disabled
               Snapshot: Post Use Sani End
        2.00
        2.00
               End Block
           Block: Pause Complete Sanitization
      0.00
        0.00
               Base: Time, ColumnSameAsMain
        0.00
               Pause: 30.00 {min}
        0.01
               End Block
      0.00
             End Block
    0.00 Block: (Storage) #Storage
             Base: SameAsMain
      0.00
      0.00
             FlowDeviation FIT PA: 20.0 {1/hour}, -20.0 {1/hour}, 300.0
      {sec}, Enabled
             FlowDeviation FIT PB: 20.0 {1/hour}, -20.0 {1/hour}, 300.0
      0.00
      {sec}, Enabled
      0.00 Air Alarm: Disabled, Enabled
      0.00 FIT_PB Totalizer Reset
      0.00 Inlet: Closed, Inlet5
      0.00 BubbleTrap: Inline
      0.00 Filter: Inline
      0.00 Column: UpFlow
      0.00 Outlet: Waste
      0.00 Flow: (300) #Storage Flowrate {cm/h}
      0.00 Set mark: Storage
      2.00
           Snapshot: Storage End
      2.00
             End Block
    0.00 Block: (Blank) #Blank
      0.00
           Base: SameAsMain
      0.00
             End Block
    0.00 Block: Return to Default
      0.00 Base: Time, Any
```

UNICORN 7.3.0 11(14)

User: AM/c304117 4/2/2025 9:45:13 AM -04:00

Method: v001 Scouting Method LB2273 CEX

Comment: This block sets skid flowpath back to default settings to

prevent nuisance valve alarms.

0.00 ManFlow: 0.0 {%}
0.05 Filter: Bypass

0.05 Column: Bypass_Both
0.10 BubbleTrap: Inline
0.10 Inlet: Closed, Closed

0.50 Outlet: Closed

0.50 End Block

0.00 Block: End of Run Delay

0.00 Base: Time, Any

Comment: This block gives system time for data transfer to server.

0.10 End Block

Scouting

Run	Included	Startup_Blocks	Prepare_Purge_Sk id Filter	<pre>Purge_Skid_Filte r Inlet 1</pre>
1	Yes	Startup_Blocks		Purge_Skid_Filte r Inlet 1
2	Yes	Blank	Blank	Blank
3	Yes	Blank	Blank	Blank
4	Yes	Blank	Blank	Blank
5	Yes	Blank	Blank	Blank
6	Yes	Blank	Blank	Blank
7	Yes	Blank	Blank	Blank
8	Yes	Blank	Blank	Blank
9	Yes	Blank	Blank	Blank
10	Yes	Blank	Blank	Blank
11	Yes	Blank	Blank	Blank
12	Yes	Blank	Blank	Blank
Run	Included	Flush_Outlet_Mai	Pause_attach_out	Pre-
Run	Included	nstreams_Equil	let_containers	Use_Sanitization
Run 1	Included Yes	nstreams_Equil Flush_Outlet_Mai	let_containers Pause_attach_out	Use_Sanitization Pre-
1	Yes	nstreams_Equil Flush_Outlet_Mai nstreams_Equil	let_containers Pause_attach_out let_containers	Use_Sanitization Pre- Use_Sanitization
1 2	Yes Yes	nstreams_Equil Flush_Outlet_Mai nstreams_Equil Blank	let_containers Pause_attach_out let_containers Blank	Use_Sanitization Pre- Use_Sanitization Blank
1 2 3	Yes Yes Yes	nstreams_Equil Flush_Outlet_Mai nstreams_Equil Blank Blank	let_containers Pause_attach_out let_containers Blank Blank	Use_Sanitization Pre- Use_Sanitization Blank Blank
1 2 3 4	Yes Yes Yes Yes	nstreams_Equil Flush_Outlet_Mai nstreams_Equil Blank Blank Blank Blank	let_containers Pause_attach_out let_containers Blank Blank Blank	Use_Sanitization Pre- Use_Sanitization Blank Blank Blank Blank
1 2 3 4 5	Yes Yes Yes Yes Yes	nstreams_Equil Flush_Outlet_Mai nstreams_Equil Blank Blank Blank Blank Blank	let_containers Pause_attach_out let_containers Blank Blank Blank Blank Blank	Use_Sanitization Pre- Use_Sanitization Blank Blank Blank Blank Blank
1 2 3 4 5 6	Yes Yes Yes Yes Yes Yes Yes	nstreams_Equil Flush_Outlet_Mai nstreams_Equil Blank Blank Blank Blank Blank Blank Blank	let_containers Pause_attach_out let_containers Blank Blank Blank Blank Blank Blank Blank	Use_Sanitization Pre- Use_Sanitization Blank Blank Blank Blank Blank Blank Blank
1 2 3 4 5 6 7	Yes Yes Yes Yes Yes	nstreams_Equil Flush_Outlet_Mai nstreams_Equil Blank Blank Blank Blank Blank	let_containers Pause_attach_out let_containers Blank Blank Blank Blank Blank	Use_Sanitization Pre- Use_Sanitization Blank Blank Blank Blank Blank
1 2 3 4 5 6 7 8	Yes	nstreams_Equil Flush_Outlet_Mai nstreams_Equil Blank	let_containers Pause_attach_out let_containers Blank	Use_Sanitization Pre- Use_Sanitization Blank
1 2 3 4 5 6 7 8 9	Yes	nstreams_Equil Flush_Outlet_Mai nstreams_Equil Blank	let_containers Pause_attach_out let_containers Blank	Use_Sanitization Pre- Use_Sanitization Blank
1 2 3 4 5 6 7 8	Yes	nstreams_Equil Flush_Outlet_Mai nstreams_Equil Blank	let_containers Pause_attach_out let_containers Blank	Use_Sanitization Pre- Use_Sanitization Blank

?

UNICORN 7.3.0 User: AM/c304117 4/2/2025 9:45:13 AM -04:00 Method: v001 Scouting Method LB2273 CEX				
Run	Included		Pause attach out	Pre-
12	Yes	nstreams_Equil Blank		Use_Sanitization Blank
12	103	DIGIIK	Diank	DIAIIK
Run	Included	Pre_Sani_Flowrat	Equil_Flowrate	Connect_Charge_t o Inlet Sample
1	Yes	300	300	Connect_Charge_t o Inlet Sample
2	Yes	300	300	Blank
3	Yes	300	300	Blank
4	Yes	300	300	Blank
5	Yes	300	300	Connect Charge t
				o Inlet Sample
6	Yes	300	300	Blank
7	Yes	300	300	Blank
8	Yes	300	300	Blank
9	Yes	300	300	Connect Charge t
				o Inlet Sample
10	Yes	300	300	 Blank
11	Yes	300	300	Blank
12	Yes	300	300	Blank
Run	Included	Charge_Flowrate		Post_Charge_Wash
1	7.7	0.00	e	_Flowrate
	Yes	200	800.00	200
2	Yes	200	800.00	200
3	Yes	200	800.00	200
4	Yes	200	800.00	200
5	Yes	200	800.00	200
6	Yes	200	800.00	200
7	Yes	200	800.00	200
8	Yes	200	800.00	200
9	Yes	200	800.00	200
10	Yes	200	800.00	200
11	Yes	200	800.00	200
12	Yes	200	800.00	200
Run	Included	Elution	Elution_Flowrate	Reduced_Elution_ Flow
1	Yes	Elution	200	150
2	Yes	Elution	200	150
3	Yes	Elution	200	150
4	Yes	Elution	200	150
5	Yes	Elution	200	150
6	Yes	Elution	200	150

?

UNICORN 7.3.0 User: AM/c304117 4/2/2025 9:45:13 AM -04:00				13 (14)
		ng Method LB2273 (
Run	Included	Elution	Elution_Flowrate	Reduced_Elution_ Flow
7	Yes	Elution	200	150
8	Yes	Elution	200	150
9	Yes	Elution	200	150
10	Yes	Elution	200	150
11	Yes	Elution	200	150
12	Yes	Blank	200	150
Run	Included	MS_Outlet	Elution_Full_Bac kside	High_Salt_Wash_F lowrate
1	Yes	Outlet1	Blank	300
2	Yes	Outlet1	Blank	300
3	Yes	Outlet1	Blank	300
4	Yes	Outlet1	Blank	300
5	Yes	Outlet2	Blank	300
6	Yes	Outlet2	Blank	300
7	Yes	Outlet2	Blank	300
8	Yes	Outlet2	Blank	300
9	Yes	Outlet3	Blank	300
10	Yes	Outlet3	Blank	300
11	Yes	Outlet3	Blank	300
12	Yes	Outlet3	Elution_Full_Bac	300
			kside	
Run	Included	Post_Use_Sani_Fl owrate	Storage	Storage_Flowrate
1	Yes	300	Blank	300
2	Yes	300	Blank	300
3	Yes	300	Blank	300
4	Yes	300	Blank	300
5	Yes	300	Blank	300
6	Yes	300	Blank	300
7	Yes	300	Blank	300
8	Yes	300	Blank	300
9	Yes	300	Blank	300
10	Yes	300	Blank	300
11	Yes	300	Blank	300
12	Yes	300	Storage	300
Run	Included	Blank		
1	Yes	Blank		
2	Yes	Blank		
3	Yes	Blank		
4	Yes	Blank		

UNICORN 7.3.0 14(14)

User: AM/c304117 4/2/2025 9:45:13 AM -04:00

Method: v001 Scouting Method LB2273 CEX

Run	Included	Blank
5	Yes	Blank
6	Yes	Blank
7	Yes	Blank
8	Yes	Blank
9	Yes	Blank
10	Yes	Blank
11	Yes	Blank
12	Yes	Blank

Method information

Signatures

Date:4/2/2025 9:44:49 AM -04:00

Description: This method is locked down for LB2273 CEX phase 1

User name:c304117
Full name:Helen Corbat

Job title: Engineer - BRD Technical Operations

Lock status: The item is locked

Questions

No. 1: Is QD00305 connected to Inlet Sample?

Question type: Mandatory
Answer type: Multiple choice

No. 2: Is QD00007 connected to Inlet 4?

Question type: Mandatory
Answer type: Multiple choice

No. 3: Is QD00008 connected to Inlet 5?

Question type: Mandatory Answer type: Multiple choice

No. 4: Is QD00346 connected to Inlet 7?

Question type: Mandatory Answer type: Multiple choice

No. 5: Is Outlet Waste directed to AWN?

Question type: Mandatory Answer type: Multiple choice

Result name and location

Result folder name: /K360PP/1 Campaigns/LB2273/Phase I/CEX

Scouting folder name: LB2273 CEX BR

Result name: LB2273 CEX BR (Type: Name) without unique identifier